




The **CEALOGIC Amniotic Perforator** is designed to puncture the amniotic membrane without any injury to the fetus or uterus. This device allows medical personnel to carefully snag the membrane without scratching or causing any damage.

Four Distinctive Features



Indication

The Recibio, Inc. **CEALOGIC Amniotic Perforator** is intended to allow medical personnel to easily and safely perforate the amniotic membrane prior to childbirth.

Additional Product Highlights

- Allows for a precise, small perforation to enable amniotic fluid extraction, if desired
- Employs a snagging not poking motion to prevent fetal injury
- Disposable and cost-effective
- Sterile

Who is Recibio, Inc.?

Recibio, Inc was founded in 2015 in Houston, Texas to explore the growing clinical application of cellular reclamation and reuse. We are committed to facilitating surgeons' ability to quickly and simply harness readily available cellular material.

Our first product offerings include **CEALOGIC Specimen Collection Kit**, **CEALOGIC Ratio Applicator Kit**, and **CEALOGIC Amniotic Perforator**.



www.recibio.us

Toll Free USA 1-844-RECIBIO

© Recibio, Inc. 2019. All rights reserved.



Made in the U.S.A. including foreign components
Assembled by **Busse Hospital Disposables**