



The **CEALOGIC Ratio Applicator Kit** from Recibio, Inc. offers surgeons a single-use option to ease simultaneous delivery of two liquids to the treatment site.

## Four Distinctive Features

### 1 Simple one-handed precision

Placement of the two syringes into a single holder/plunger creates a uniform mixture of the two fluids that can then be precisely sprayed onto the treatment site.

### 2 Finely micronized spray


Spray tip with Dual-Lock Adapter produces an exceptionally small droplet size (20 – 230 microns) to form a highly concentrated mist.

### 3 Meets elite performance standards

Dispenses 13mL in total deliverable volume with optimized flow rate, droplet size, droplet velocity, and coverage.

### 4 Steady stream sprayer tip for accuracy

Generates a narrow, concentrated spray plume to ensure more accurate coverage of the treatment site.



# CEALOGIC™ Ratio Applicator Kit


All components are shipped sterile in package.

## Indication

The components of the Recibio, Inc. **CEALOGIC Ratio Applicator Kit** are intended to assist the user in simultaneously delivering two liquids to the treatment site.

## Additional Product Highlights

- Kit components are readily assembled to allow mixing of two fluids at an 12:1 ratio
- Each kit comes with two standard syringes (one 12 cc, one 1 cc) that are easily placed into the dispenser
- Syringes are transparent for unobstructed view of fluids and calibrated to ensure precise measurement
- Each kit comes with cup and label to send for analysis, if desired
- All components are latex free



## Each kit contains the following components:

- One Sprayer Tip
- One Sprayer Dispenser
- One Sprayer Adapter
- One Specimen Cup
- One Specimen Label
- One Plastic Cup with Foam Liner
- One 12cc Luer Lock Syringe
- One 1cc Luer Lock Syringe
- One Assembly Diagram/information Booklet
- One Product Sticker Sheet

## Who is Recibio, Inc.?

Recibio, Inc was founded in 2015 in Houston, Texas to explore the growing clinical application of cellular reclamation and reuse. We are committed to facilitating surgeons' ability to quickly and simply harness readily available cellular material.

Our first product offerings include **CEALOGIC Specimen Collection Kit**, **CEALOGIC Ratio Applicator Kit**, and **CEALOGIC Amniotic Perforator**.



[www.recibio.us](http://www.recibio.us)

Toll Free USA 1-844-RECIBIO

© Recibio, Inc. 2019. All rights reserved.



[www.cealogic.com](http://www.cealogic.com)

Made in the U.S.A. including foreign components  
Assembled by **Busse Hospital Disposables**