



The **CEALOGIC Specimen Collection Kit** from Recibio, Inc. has a unique, patent-pending design to ease both specimen collection and removal. The potential uses for this specific collection container are wide-ranging.

Four Distinctive Features

1 Patent-pending, angled, Luer lock™ port

The port to draw fluid from the specimen container is uniquely set at a 45-degree angle to ease attachment of the syringe.

2 Innovative specimen removal port

The blue-capped removal port allows the specimen to be drawn directly from the collection container to preserve sterility during transfer.

3 Unobstructed specimen evaluation

The container is transparent to ensure clear specimen visualization and assessment.

4 Large container capacity

The 80cc calibrated container allows for maximum specimen collection and ensures precise measurement.



CEALOGIC™ Specimen Collection Kit Package

All components are shipped sterile in package.

Indication

The Recibio, Inc. **CEALOGIC Specimen Collection Kit** is intended for use of normal suction and specimen collection.

Additional Product Highlights

- All-in-one convenience packaging eases specimen collection and ensures sterility
- Each kit comes with cap and labels to send to analysis, if desired
- Narrow, round-tip Yankauer suction promotes safety during collection
- 10-foot tubing and 5-in-1 connector allows for connection to any vacuum
- All components are latex free



Each kit contains the following components:

- One Specimen Collection Container
- One 10'-1/4" Suction Tubing with Male Adapter
- One Yankauer Suction Tip
- One Specimen Cap
- One Specimen Label
- One Product Sticker Sheet
- One Assembly Diagram/Information Booklet
- One 5-in-1 Connector

Who is Recibio, Inc.?

Recibio, Inc was founded in 2015 in Houston, Texas to explore the growing clinical application of cellular reclamation and reuse. We are committed to facilitating surgeons' ability to quickly and simply harness readily available cellular material.

Our first product offerings include **CEALOGIC Specimen Collection Kit**, **CEALOGIC Ratio Applicator Kit**, and **CEALOGIC Amniotic Perforator**.



www.recibio.us

Toll Free USA 1-844-RECIBIO

© Recibio, Inc. 2019. All rights reserved.



www.cealogic.com

Made in the U.S.A. including foreign components
Assembled by **Busse Hospital Disposables**